



## Static Human Blood Clot Tests (Midline Catheters)

**Comparing the Hydromer 2018-20M and Hydromer EF50L coatings to an uncoated Midline catheter.**

**Presented by Martin von Dyck  
(Executive VP-Hydromer)  
[www.Hydromer.com](http://www.Hydromer.com)**



## “STATIC HUMAN BLOOD CLOT TEST”


- **OBJECTIVE**: The objective of this test is to compare the gross thrombo-resistivity of a Midline (made using using 95A Tecothane®) catheters when soaking in re-calcified citrated human blood.
- After 5 hours soaking in the blood at 37°C the catheters will be gently rinsed in PBS and photographed.
- This test will provide a gross indication of how the respective catheters and coatings should perform in the more detailed testing to follow.

## “STATIC HUMAN BLOOD CLOT TEST” Protocol

- Midline catheters were coated with Hydromer 2018-20M and Hydromer EF50L as per Hydromer’s coating procedures.
- A Crystal Violet stain test will confirmed the presence and uniformity of the Hydromer coatings.
- *All* catheter samples (uncoated or coated) will be incubated in re-calcified citrated human whole blood for 5 hours @ 37<sup>0</sup>C and then gently rinsed with PBS and photographed.
- *Blood Source:* Human Unspun Whole Blood with Sodium Citrate was drawn from a male on July 31, 2017. (4 days prior to testing.)
- *Blood Prep.:* Blood recalcification was performed by Hydromer R&D using 200ml Citrated Whole Human Blood and adding 0.92ml of CaCl<sub>2</sub> resulting in a 4 milliMolar final concentration. Complete clot formation was verified at 45min.

# Hydromer coating “Crystal Violet stain test” to confirm coating uniformity.


Crystal Violet stain confirmed the presence of the coating distal from the suture junction.



2018-20M

EF50L

Crystal Violet stain confirmed the presence of the coating in the lumens.



2018-20M


EF50L



# Static Blood Clot Test set-up



Sample tested  
#1 Uncoated Control  
#2 Hydromer 2018-20M  
#3 Hydromer EF50L



# Static Blood Clot Test set-up



**Positive control verifying that the blood clotted before the 5 hour test limit. Actual clotting occurred at 45 minutes.**

# Static Blood Clot Test

## Uncoated Midline catheter


Catheters incubated in re-calcified citrated human whole blood.  
(For 5 hours at 37°C)




Uncoated Catheter was completely covered by clotted blood. Distal end of blood clot contracted during removal from glass test tube exposing the distal end of the Uncoated catheter.

# Static Blood Clot Test Uncoated Midline catheter

Catheters incubated in re-calcified citrated human whole blood.  
(For 5 hours at 37°C)



**Clot Sheath was physically removed from Uncoated Midline and flushed.**





# Static Blood Clot Test

## Hydromer 2018-20M coated Midline catheter

Catheters incubated in re-calcified citrated human whole blood.  
(For 5 hours at 37°C)



The blood clot slipped off of the Hydromer coated catheter showing no residual clot bonds on the catheter.



Lumens are clear!  
Exterior is clear!




# Static Blood Clot Test

## Hydromer EF50L coated Midline catheter

Catheters incubated in re-calcified citrated human whole blood.



(For 5 hours at 37°C)




The blood clot slipped off of the Hydromer coated catheter showing no residual clot bonds on the catheter.



Lumens are clear!  
Exterior is clear!



**Catheters incubated in re-calcified citrated human whole blood.  
(For 5 hours at 37°C)**



## Conclusion

Based on this in-vitro laboratory Static Blood Clot Test, both of the Hydromer coating candidates, EF50L and 2018-20M, resisted blood clot adhesion on polyurethane Midline catheters.